

## HOW TO REMPLACE THE GERMIBLOC CLEANING AND DESINFECTANT PRODUCT ONCE THE SIGNAL BLEEPS 3 TIMES AT THE END OF IT'S CYCLE.

- 1. Remove the cover with the key and turn the switch off.
- 2. Unscrew the top of the mixing product.
- 3. Replace the empty bottle by carefully screwing the new container on the bottle top.
- 4. Close the mixer (tighten it enough to lock it in perfectly).
- 5. Turn on the switch (3 bleeps) then place the cover back on.











Box with 4 Germibloc™ products = SUP1090

## HOW TO REPLACE THE BLUE CLEANING PRODUCT FOR THE SUP1080 SPRAY OPTION



Firstly, pour 250 ml of Biosurf



Then, top it up with water. Around 2,5 L



A blue cleaning product is sprayed when the user arrives for 4"

Do the same transaction once the spray in the toilet pan is no longer blue.

Box with 18 blue cleaning products = SUP1085